

# Skin Blends TOUCH100

## Electrocoagulation System

**For the quick and  
easy removal of:**

- **Distended Capillaries**
- **Spider Veins**
- **Cherry Angiomas**
- **Skin Tags**
- **Vascular Blemishes**



**One Foot Pedal  
Four Functions:**

- **Radio Frequency (RF)**
- **Galvanic**
- **Phoresis**
- **Air – Cools/Numbs**



**Before**



**After**

**You Don't Need A Laser!**

**Contact Skin Blends Today to Order  
1-877-754-6253 (1-UPS-SKINBLENDS)**


**[www.SkinBlends.com](http://www.SkinBlends.com)**

### **Cherry Angioma**

**Before**



**After**



### **Telangiectasias**




### **Vascular Blemish**




### **Skin Tags**

**Before**



**After**



### **Skin Blends TOUCH100™ is Non-invasive**

TOUCH100™ is an inexpensive and indispensable tool which allow skin care professionals to quickly and easily remove common vascular blemishes and skin tags. It yields superior results when compared to other electronic devices used by medical personnel that commonly cause scarring, hypo-pigmentation, and hyper-pigmentation. When compared to laser treatment, the TOUCH100™ produces similar results without the pain, bruising, or high treatment costs.

### **Why TOUCH100™ Works Better**

This quick and simple treatment will not cause scarring or hyperpigmentation. TOUCH100™ uses very low levels of Radio Frequency (RF) and DC galvanic current simultaneously to electrocoagulate the vessel, without desiccating surrounding tissue. This treatment is self-limiting and vessel specific, resulting in minimal skin reaction. Low treatment energy does not adversely affect the dermis to cause scarring or any other major reaction.

### **How Long Does the Treatment Take?**

Most blemishes are removed in seconds. Depending on the number of capillaries, most conditions can be successfully and permanently corrected in one or two 10–15 minute sessions.

### **How is the Non-Invasive Procedure Done?**

The surface of the skin is cleansed and prepped. A fine tipped sterile probe is placed on the surface of the skin. There is NO pressure or penetration of the skin during treatment. Low-level RF and DC (direct current) are

passed through the probe to the vessel. The tiny vessel is clotted or coagulated and blood stops flowing. With proper post treatment care the vessel is naturally absorbed into the surrounding tissue. Typically within three to four weeks the vessel disappears. Most skin tags fall off within days of being treated.

### **Includes Phoresis Mode Which:**

- Calms skin and reduces swelling;
- Increases penetration of “ampule” and other products;
- Facilitates disencrustation treatments.

### **Includes Air Desensitization Mode**

Directing a stream of blowing air at the treatment site fools nerves in the skin, creating an anesthetic effect that minimizes any patient discomfort.

### **High Technology at Low Cost**

Treatment procedures are easy to perform, and generate on average between \$75 and \$150 per treatment. The TOUCH100™ pays for itself quickly, creating new business and satisfied clients.

Contact Skin Blends Today to Order:  
1-877-754-6253

FREE training included  
with purchase of  
system!

**NOTE: Each state has a regulatory authority responsible for governing the use of electrical probe devices by a nonmedical therapist. You should inquire of your state regulatory agency as to any determination regarding treatment with this type device.**